

Below are the instructions on how to take care of your Silkscin Pads:

PRIOR TO USE

Prior to use, skin should be absolutely free of any products i.e. moisturizers, soaps, oils, soap residue from the shower / bath, spf's etc. These can all affect the adhesive properties of the Pad and either reduce its longevity or render it ineffective.

WASHING

When washing the pads, we recommend the use of either Sunlight Laundry soap (as it's sometimes known) as it contains no moisturizers. Should you use soap which does contain moisturizers (as most do), this will compromise the integrity of the silicone therefore reducing the use you will get from your pad.

The Sunlight Laundry soap is easily obtainable from the laundry section of the supermarket. No costly soaps or serums are required to wash the pads - laundry soap lasts a very long time, is very economical, and washes the pads perfectly

HOW OFTEN TO WASH

The Pads should be washed in warm soapy water after the first or second use and as required thereafter (approximately weekly) i.e. when there is a noticeable build-up of skin cells and atmospheric ingress on the pads surface.

When the weather is warmer, the skin tends to sweat more hence the potential need for more frequent washing (every 5-7 days). Once washed, place the pads sticky side up on the plastic plate until dry then simply turn it over so it adheres to the plate and is ready for use.